

KRUUSE Rehab Knee Protection

Size Chart

Size	Cat. No	Suitable for
XXX-Small, left	279944	Chihuahua, Japanese Chin, Maltese, Pomeranian, Japanese Terrier, Yorkshire Terrier, Papillon
XXX-Small, right	279945	
XX-Small, left	279942	Bichon Frisé, Coton de Tulear, Shih Tzu, Tibetan Spaniel, English Toy Spaniel, Yorkshire terrier
XX-Small, right	279943	
X-Small, left	279849	Sheltie, Cocker Spaniel, Fox Terrier, French Bulldogs, Miniature Schnauzer, Cavalier king Charles Spaniel
X-Small, right	279850	
Small, left	279851	Kelpie, Whippet, Bulldog, Beagle, Dachshund, Pembroke Welsh Corgi, Boston Terrier
Small, right	279852	
Medium, left	279853	Labrador, Border Collie, Springer Spaniel, Labrador Retriever, Golden Retriever, Poodle, Boxer, Pointer (German Shorthaired), Doberman, Australian Shepherd
Medium, right	279854	
Large, left	279855	German Shepard, Boxer, Rhodesian Ridgeback, Siberian Husky, Belgian Malinois, Alaskan Malamute
Large, right	279856	
X-Large, left	279857	Rottweiler, Bullmastiff, Bloodhound, Dogue de Bordeaux, Cane Corso, Akita
X-Large, right	279858	
XX-Large, left	279885	Bernese Mountain Dog, Great Dane, Leonberger, Irish Wolfhound, Mastiff, Saint Bernard
XX-Large, right	279886	



	Size chart KRUUSE Rehab Knee Brace in CM	XXXS	XXS	XS	S	M	L	XL	XXL
A	Total length, along fold edge	7	12	20.5	25	26.5	26.5	33	33
B	Width, top edge	6.5	7.5	11.5	13.5	14.5	15.5	17.5	18.5
C	Width, (0.5 cm)* 1 cm above knee dart	5*	6*	8	10	10.5	11.5	12	13
D	Width, bottom hem	4.5	5	5	5.5	6	6.5	7	8
E	Width, widest point	7	9	15	16.5	17.5	18.5	20.5	21.5
	Strap full length	50	55	60	61	90	91	96	96.5

